

GloMax® Instruments

Perfect Partnership of Assays and Instruments

Use the instruments designed for your Promega assays

- → Select your Promega assay from the pre-loaded protocols
- → Start collecting the results.
 - No Guessing
 - No Gain Settings

What could be easier!

Get FREE Assays with Your GloMax® Instrument Purchase

Choose:

- 3 Small Assay Kits
- OR -
- 1 Large Assay Kit

See the back of this flyer for more information!



GloMax® Discover

Luminescence ● Fluorescence ●
Absorbance ● FRET ● BRET ●
6- to 384-well



GloMax® Explorer

Luminescence ● Fluorescence ● with or without Visible Absorbance ● 6– to 384-well



GloMax® Navigator
Luminescence • 96-well



GloMax® 20/20

Luminescence • Fluorescence • Single Tube

GloMax® Multi Jr.

Luminescence ● Fluorescence ● Absorbance ● Single Tube



Buy a GloMax[®] Instrument & Get up to 3 assays FREE

Choose from assays				Ohaaa		
like these:		3 Small Kits		Choos –OR–	_	
	Tomas			-UK-		
Reporter Assay	Туре	Cat.#	Size		Cat.#	Size
Nano-Glo® Dual-Luciferase® Reporter Assay System	CBA	N1610	10ml		N1620	100ml
Nano-Glo® Luciferase Assay System	CBA	N1110	10ml		N1120	10 x 10ml
ONE-Glo™ EX Luciferase Assay System	CBA	E8110	10ml		E8120	100ml
Cell Viability						
RealTime-Glo™ MT Cell Viability Assay	CBA	G9711	100 rxn		G9713	1,000 rxn
CellTiter-Glo® 2.0 Assay	CBA	G9241	10ml		G9242	100ml
CellTiter-Glo® 3D Cell Viability Assay	CBA	G9681	10ml		G9683	100ml
Apoptosis						
Caspase-Glo® 3/7 Assay	CBA	G8091	10ml		G8092	100ml
Oxidative Stress						
GSH/GSSG-Glo™ Assay	CBA	V6611	10ml		V6612	50ml
ROS-Glo™ H ₂ O ₂ Assay	CBA	G8820	10ml		G8821	50ml
Metabolism						
NAD/NADH-Glo™ Assay	CBA	G9071	10ml		G9072	50ml
NADP/NADPH-Glo™ Assay	CBA	G9081	10ml		G8092	50ml
Kinase Assays						
ADP-Glo™ Kinase Assay	ВСН				V9101	1000 assays
Kinase Enzyme System + ADP-Glo™ Kinase Assay	ВСН				>170	Choices
Multiplexed Fluorescent/Luminescence Assay						
ApoTox-Glo™ Triplex Assay	CBA				G6320	10ml
Mitochondrial ToxGlo™ Assay	СВА	G8000	10ml			
Fluorescent Assay (multiplex capable)						
CellTox™ Green Cytotoxicity Assay	CBA	G8741	10ml		G8742	50ml
CellTiter-Fluor™ Cell Viability Assay	CBA	G6080	10ml		G6081	5 x 10ml

CBA = Cell-Based Assay; BCH = In vitro Biochemical Assay

Dozens more assays to choose from... Ask us!

To receive your FREE assays, provide your Promega Representative with the catalog number(s) and your shipping address. Please allow 2--3 weeks for delivery.

ADP-Glo, ApoTox-Glo, CellTiter-Fluor, CellTox, GSH/GSSG-Glo, NAD/NADH-Glo, NADP/NADP-Glo, ONE-Glo, ROS-Glo and ToxGlo are trademarks; Caspase-Glo, Celltiter-Glo, Dual-Glo, GloMax, and Nano-Glo are registered trademarks of Promega Corporation. Products may be covered by pending or issued patents or may have certain limitation. Please visit our website for more information.

