

Labnet Prism™ Mini Centrifuge



- Marca Labnet
- Para 8 microtubos de 1.5 y 2.0 mL y 4 strips de 8 x 0.2 mL
- Rotores intercambiables tipo *quick release*
- 6000 rpm
- Frenado rápido
- Libre de mantenimiento

Especificaciones:

Capacity	8 x 1.5/2.0 mL tube rotor & 4 x 0,2 mL strips
Maximum Speed	6000 RPM
Maximum RCF	2000 x G
Dimensions	5.2 x 5.2 x 4.4 inch / 13.2 x 13.2 x 11.2 cm
Weight	1.2lbs / .55 kg
Electrical	120V or 230V, 50/50Hz

Accesorios:

SKU	TITLE
C1206	0.4/0.25 ml Adapters, Pack of 6
C1205	0.5/0.6 ml tube Adapters, Pack of 6
C1222	0.2 ml Thermal Cycling Tube Adapters, Pack of 6
C1206-8	Adapters for 0.4 mL tubes, pk 8
C1205-8	Adapters for 0.5 mL tubes, pk 8
C1222-8	Adapters for 0.2 mL tubes, pk 8