

PIPETTE	STANDARD TIPS						FILTER BARRIER TIPS						PIPETTE
	10µl Short T-10	20µl Extended T-10XL	200µl T-200	300µl T-300	1000µl T-1000	1250µl Extended T-1250	10µl Short FT-10	10µl Extended FT-10XL	20µl FT-20	100µl FT-100	200µl FT-200	1000µl FT-1000	
<b>SARTORIUS BIOHIT</b>													
Picus e-Line Single, 8 & 12 channel	0.1-2.5µl	0.1-10µl	0.2-10µl	5-120µl	50-1000µl		0.1-2.5µl	0.1-10µl					
ProLine ProLine Plus m-Line Single channel	0.1-3µl	0.2-10µl	0.5-10µl	5-50µl	5-100µl	10-250µl	20-200µl	50-1000µl					
ProLine ProLine Plus m-Line 8 and 12 channel	0.5-10µl	5-100µl	25-250µl (ProLine)	50-250µl (ProLine)	30-300µl (mLine)								
<b>BRANDTECH</b>													
Transferpette Single channel	0.5-10µl	2-20µl	10-100µl	20-200µl	25-250µl	100-1000µl							
Transferpette 8 and 12 channel	0.5-10µl	2-20µl	2.5-25µl	10-100µl	20-200µl								
<b>CAPP</b>													
AERO Single, 8 and 12 channel	0.2-2µl	0.5-10µl	10-100µl	25-200µl	50-300µl	100-1000µl							
<b>EPPENDORF</b>													
Research Single, Multi channel	0.1-2.5µl	0.5-10µl	2-20µl	10-100µl	20-200µl	100-1000µl	20-300µl (m/ch)	30-300µl (m/ch)					
Reference Single, Multi channel	0.1-2.5µl	0.5-10µl	2-20µl (UltraMicro)	2-20µl	10-100µl	20-200µl	50-200µl	100-1000µl					
Electronic Single, Multi channel	0.5-10µl	5-100µl	2-200µl	50-1000µl									
<b>THERMO FINNPIPETTE</b>													
0.2-2µl (Micro) Finnpipette Single channel													
1-10µl (Micro)													
0.5-10µl													
10-100µl													
20-200µl													
30-300µl													
100-1000µl													
1-10µl Finn Multi channel													
50-300µl													
0.5-5µl (Micro)													
1-10µl (Micro) Focus Single channel													
3-30µl													
5-50µl													
10-100µl													
30-300µl													
100-1000µl													
1-10µl Electronic Multi channel													
30-300µl													
100-1200µl													
<b>GILSON</b>													
0.2-2µl, P2/P2N Pipetman Classic/Neo Single channel													
1-10µl, P10/P10N													
2-20µl, P20/P20N													
20-100µl, P100/P100N													
50-200µl, P200/P200N													
200-1000µl, P1000/P1000N													
<b>HTL LABMATE, DISCOVERY*</b>													
0.2-2µl Single channel													
0.5-10µl													
2-20µl													
10-100µl													
20-200µl													
100-1000µl													
1-10µl Multi channel													
5-50µl													
20-200µl													
50-300µl													
<b>METTLER TOLEDO/RAININ</b>													
0.1-2µl Pipet-Lite/XLS (standard shafts) Single, Multi channel													
0.2-10µl													
2-20µl													
10-100µl													
20-200µl													
100-1000µl													
<b>NICHIRYO</b>													
0.5-10µl Nichipet EX / Oxford Single channel													
2-20µl													
10-100µl													
40-200µl													
100-1000µl													
50-300µl, m/ch													
<b>SOCOREX</b>													
0.1-2µl Accura Single channel													
0.5-10µl, single + m/ch													
2-20µl													
10-100µl, single + m/ch													
20-200µl, single + m/ch													
100-1000µl													

This table is intended to be a general guide only. Pipette manufacturers occasionally change nose cone/shaft design and so tip fit can not be guaranteed. \* includes brands: Axygen AxyPet, Labnet Biopipette, Genemate, Cleaver Scientific Omnipipette, VWR Signature