

Novas Bio Compatibility Chart : PCR Plates

Novas Part No.	Instrument Manufacturer																	
	Agilent®		Analytik Jena®				Biorad®							Eppendorf®	Techne®			
	AriaMx	Robocycler Gradient	qTOWER ³	Alpha SC	Flexcycler/2	Biometra TOne/TAAdvanced/Gradient/Pro (96 Series)	Biometra TRobot II	CFX Connect	CFX 96/96 Touch/96 Touch Deep Well	Chromo4	C1000/C1000 Touch	S1000	MJ Research DNA Engine Opticon/Opticon 2	MJ Research DNA Engine Tetrad/Tetrad 2	MJ Research DNA Engine Dyad/Dyad Disciple	MJ Research PTC 100/200/220/221/225/240	Mastercycler ep realplex	Quantica
HP96-LPFS-WW	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•
HP96-LPFS-BKC	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•

LEGEND	
Real Time Instruments (qPCR)	
Thermocyclers (PCR)	